

FDPH <hybrid>

Fully Digital Hybrid Flatpanel Solutions

for Chiropractic, Veterinarian, & Private Practice offices

Taking digital radiography to the next level with our 17x17 Hybrid Flat panel technology is 20/20 Imaging's cassette-sized FPH. Enjoy 3 unique operating configurations, including full wireless, wired power connection, and wired (fixed) power and data connections. Providing the most in versatility with exceptional image quality, and paired with advanced software features and functionality. The 17x17 hybrid panel uses enhanced Gadolinium technology and has many unique features such as device-drop notification emails.

The FPH is ideal for use with both new x-ray rooms as well as a retrofit solution. The versatile panel utilizes the Opal-RAD Acquire software for image capture, manipulation, and sending. 20/20's acquire software offers useful tools to enhance the image acquisition process, thus keeping your office in the forefront of x-ray technology.

SPECIFICATIONS

- Warranty: 5 years
- Imaging Area: 17" x 17"
- Line pair per mm: 3.6
- Wired/Wireless connectivity
- Cassette sized
- Battery Life: >3.5hrs
- 9.5 lbs
- Pixel Pitch: 139um
- Off-center Imaging
- Dim: 18.11"x18.11"x0.6"
- 460x460x15.1mm
- A/D conversion: 14-bit

3 unique operating configurations:
data+power, power-only, wireless

OPAL-CHIRO



OPAL-VET



Customized Practice-Specific Software

20/20 Imaging integrates the specialized Opal Chiropractic tool-set for the Chiropractic discipline offering advanced measurement and alignment tools, and Veterinary specific animal configuration for the Vet specialty. Every New or Retro-fit 2020 system comes complete with a set of advanced PACS tools and functionalities for your office.

Cassette Sized. Fits any existing system

The FPH is cassette sized, giving you the ability to add to any x-ray unit. 20/20 Imaging also offers new x-ray units for practices without existing x-ray, or those who are simply looking to upgrade their older x-ray equipment.



20/20's Straight Arm DR (left) offers an advanced digital x-ray solution in an efficient and compact design. It allows for all imaging views that are required in radiology including wheelchair and table work for maximum patient flexibility. Retro-fit solutions for existing x-ray units also available (right).



20/20's New Digital X-Ray System V-FPH (left) comes complete with Acquisition Interface touch-screen built in. The V-FPH also features a high-frequency generator providing superior image quality. Retro-fit solutions for existing x-ray units also available (right).



471 Jennings Dr. Ste A Lake in the Hills, IL 60156

p: 866-734-6234 • 2020imaging.net/contact