

FPW <WIRELESS>

Fully Digital Wireless Flatpanel Solutions

for Chiropractic, Veterinarian, & Private Practice offices

Taking digital radiography to the next level with our Wireless Flat panel technology is 20/20 Imaging's cassette-sized FPW. Providing exceptional image quality, paired with advanced software features and functionality. The 14x17 panel uses enhanced Gadolinium technology and has many unique features such as device-drop notification emails.

The FPW is ideal for use with both new x-ray rooms as well as a retrofit solution. The versatile panel utilizes the Opal-RAD Acquire software for image capture, manipulation, and sending. 20/20's acquire software offers useful tools to enhance the image acquisition process, thus keeping your office in the forefront of x-ray technology.

SPECIFICATIONS

- Warranty: 3 years
- Imaging Area: 14" x 17"
- Line pair per mm: 3.3
- Wireless connectivity
- Cassette sized
- Battery Life: >3.5hrs
- 7.93 lbs
- Pixel Matrix: 2304x2800
- Off-center Imaging
- Dimensions: 15.12"x18.11"x0.6"
- A/D conversion: 14-bit

Cassette Sized. Fits any existing system

The FPW is cassette sized, giving you the ability to add to any x-ray unit. 20/20 Imaging also offers new x-ray units for practices without existing x-ray, or those who are simply looking to upgrade their older x-ray equipment.

OPAL-CHIRO



OPAL-VET



Customized Practice-Specific Software

20/20 Imaging integrates the specialized Viztek chiropractic tool-set for the Chiropractic discipline offering advanced measurement and alignment tools, and Veterinary specific animal configuration for the Vet specialty. Every New or Retro-fit 20/20 system comes complete with a set of advanced PACS tools and functionalities for your office.



20/20's Straight Arm DR (left) offers an advanced digital x-ray solution in an efficient and compact design. It allows for all imaging views that are required in radiology including wheelchair and table work for maximum patient flexibility. Retro-fit solutions for existing x-ray units also available (right).



20/20's New Digital X-Ray System V-FPW (left) comes complete with Acquisition Interface touch-screen built in. The V-FPW also features a high-frequency generator providing superior image quality. Retro-fit solutions for existing x-ray units also available (right).



a division of Konica Minolta Medical Imaging

KONICA MINOLTA

471 Jennings Dr. Ste A Lake in the Hills, IL 60156

p: 866-734-6234 • 2020imaging.net/contact