



DIGITAL X-RAY PODIATRY SOLUTIONS

Revolutionary Imaging Systems



2020imaging.net/podiatry

Wireless and Tethered DR

MOMENTUM DR 10X12 WIRELESS PANEL



MFG BY KONICA MINOLTA

- Wireless Interface
- Featherlight 3.30 lbs
- Cesium Scintillator
- 100 micron | 5 lp/mm
- Lithium Ion Capacitor (7-yr warranty)
- Cassette Sized
- 16-bit Images
- Rapid Charge

NEW UNIVERSAL DR 10X13 WIRELESS PANEL



UNIVERSAL DR

- Wireless Interface
- Lightweight 5 lbs
- Cesium Scintillator
- 100 micron | 5 lp/mm
- Lithium Ion Battery **Extended Battery Life**
- Active Image Area: 10"x13"
- 16-bit Images

UNIVERSAL DR 1012V TETHERED PANEL



UNIVERSAL DR

- Tethered Interface
- Lightweight 5.5 lbs
- Cesium Scintillator
- 125 micron | 4 lp/mm
- Detachable & Durable Cable
- Cassette Sized
- 16-bit Images

UNIVERSAL DR 14X17 WIRELESS PANEL



UNIVERSAL DR

- Wireless Interface
- Lightweight 7.3 lbs
- Cesium Scintillator
- 150 micron | 3.4 lp/mm
- Lithium Ion Battery
- Cassette Sized
- 16-bit Images

Optional high resolution panel available

New and Retrofit X-Ray

NEW PXS710D X-RAY

- High Frequency Output*
- Dose Efficient/Ext Tube Life
- Superb Image Quality
- Bi-directional Tube Movement
- Integrated Cable Management
- Enhanced Patient Support
- UL Certified 450lbs weight limit
- Qualifies for ADA Tax Credit
- UL & CUL Listed
- Health Canada Approved

New X-Cel X-Ray models also available

RETROFIT ANY X-RAY

Wireless & Tethered DR Compatible
All options qualify for ADA Tax Credit



Example of older X-cel MB-700 modified with new 3.5'h SD Base and 1012V panel.

*Compared to standard frequency (single phase) x-ray.

Hospital-Grade Software

ACQUISITION INTERFACE



- Under 6 second ⌚ Processing
- Preset X-ray Exams
- Anatomical Program (APR)
- Advanced Image Processing
- Touchscreen Enabled
- Background Submission
- Adult & Pediatric Settings
- Fully Customizable
- Image Enhancement Tools
- User Friendly Interface

ADVANCED OPAL-RAD VIEWER



Multi-funtion Studylist ⬇️ Advanced Viewer ➡️

- EHR Integration Ready**
- Burn/Import CDs
- Custom Search/Query
- Customized User Settings
- Podiatric Toolset Available†
- Fast, User Friendly Interface
- Image Fine Tuning
- Exam Room Viewing



MOBILE VIEWING WITH M-OPAL†

iPad®, iPhone®, Android™, Blackberry®, any mobile or smart device!



- Lightweight & Fast
- Patient List & Viewer
- Simple Annotations Enabled
- Tap, Drag, and Pinch
- Patient Filter & Search



CLOUD-BACKUP OPTIONAL

- Live Synchronization
- Fully Monitored
- 448-bit Encryption (traffic)
- Fully Customizable
- Enterprise Grade
- Recover Deleted Files



Advanced
Superior
Efficient
Reliable
Affordable

2020imaging.net/podiatry

