

**md buyline**      **Konica Minolta**  
**#1 in DR Customer Satisfaction**

Market Intelligence Briefing™ Q1 2017, Radiographic Unit, Digital and Conventional



panel docked in charging station

20/20 Imaging introduces the newest member of the DR family for Podiatrists. The 10x12 inch wireless flat panel detector is lightweight and easy to maneuver, and is a true cassette sized DR Panel with a Cesium scintillator; providing the highest image quality possible.

The DRMW is ideal for use with both new x-ray as well as a retrofit solutions. The versatile panel utilizes the Opal-RAD Acquire software for image capture, manipulation, and sending. 2020's acquisition software offers useful tools to enhance the image acquisition process, thus keeping your office in the forefront of x-ray technology.

### SPECIFICATIONS

- Imaging Area: 10" x 12"
- Scintillator: CsI (Cesium)
- Line Pair: 3.4lp/mm
- Wireless connectivity
- Cassette sized
- Rapid Charge enabled
- Battery Life: 150 images/4hrs
- 16-bit AD Conversion
- 3.75 lbs
- Dimensions: 11.09"x13.11"x0.626"

[2020imaging.net/momentum](http://2020imaging.net/momentum)



### Customized Practice-Specific Software\*

20/20 Imaging integrates the specialized Opal Podiatry tool-set for the Podiatric discipline offering advanced measurement and alignment tools. Every New or Retro-fit 2020 system comes complete with a set of advanced PACS tools and functionalities for your office.

\*Available for an additional cost. Contact sales for details.

