

md buyline

**Konica Minolta
#1 in DR Customer Satisfaction**

Market Intelligence Briefing™ Q2 2017, Radiographic Unit, Digital and Conventional



panel shown docked in charging station

Manufactured by Konica Minolta, the MOMENTUM-DR 10x12 inch wireless flat panel detector is lightweight and easy to maneuver, and is a true cassette sized DR Panel with a Cesium scintillator; providing the highest image quality possible. The Momentum is ideal for use with both new x-ray as well as a retrofit solutions.

Ranked #1 in DR Customer Satisfaction for the last 7 quarters; OEM's now choose Konica panels as their "go-to" panel over their own panels! Konica panels have become an industry standard of durability and image quality.

Saving your office time, the Momentum utilizes built-in **Smart Capacitors** instead of batteries, which allows for over a million charge cycles, ultra fast charging, and leaves you worry-free!

2020imaging.net/momentum

SPECIFICATIONS

- Imaging Area: 10" x 12"
- Scintillator: CsI (Cesium)
- Line Pair: 3.4lp/mm
- Wireless connectivity
- Cassette sized
- Capacitor Charging Technology
- Rapid Charge enabled
- 16-bit AD Conversion
- 3.75 lbs
- 3yr Warranty (7-yr capacitors)



Customized Practice-Specific Software*

20/20 Imaging integrates the specialized Opal Podiatry tool-set for the Podiatric discipline offering advanced measurement and alignment tools. Every New or Retro-fit 2020 system comes complete with a set of advanced PACS tools and functionalities for your office.

*Available for an additional cost. Contact sales for details.

