


# pDR HG+

## 715A-BD - 10x12 HG+ Complete Podiatry X-Ray System



### BD-715 X-RAY Features

- Versatile sensor positioning (improve patient throughput)
- Bi-Lateral view positioning (eliminates unsafe patient maneuvering)
- 3.5in. low base (improved height for patient use)
- Patient support rail system (allows access from side or front)
- 10.25"x12.75" active image matrix

### Qualifies for ADA Tax Credit

- X-ray designed and manufactured by X-Cel X-Ray exclusively for 2020 Imaging



### IMAGEVIEWER

### PATIENT STUDYLIST

USER FRIENDLY VIEWER  
DPM SPECIFIC TOOLS  
DICOM COMPLIANT IMAGES

MULTI FUNCTIONAL  
*burn & import cds  
compare images (pre/post op)  
custom search/query  
customize user privileges*



### pDR FEATURES

**NEW 2 YEAR WARRANTY**  
*upgraded from the previous 1yr warranty!*

### AUTO PROCESSING

- Auto Cropping •Auto Contrast (W/L)

### EASY MANEUVERING

- Sleek, Compact Design
- Redesigned Cable Management •Multi-Handle Configuration

### SUPERB X-RAY DETAIL

- Superior Image Quality •High Sensitivity •High Resolution Output
- 3.4 lp/mm •Flat Panel Technology

### ACQUIRE INTERFACE



**AUTOMATIC CONTRAST**  
*UPON PROCESSING*

**IMAGE AUTO SHUTTER**  
*automatically crops each view*

**BACKGROUND TRASMIT**  
*automatically sends images during study*

**X-PIPE PROCESSING TECHNOLOGY**  
*sharper images, better diagnosis*



471 Jennings Dr. Ste A Lake in the Hills, IL 60156  
Phone: 866-734-6234 • Fax: 919.772.2810  
Email: info@2020imaging.net • www.2020imaging.net